



# Procarta™ Cytokine Assays


## High-throughput Multiplex Immunoassays

### Procarta Cytokine Assay Kits


Procarta Cytokine assays simultaneously quantitate 1-28 cytokines from diverse matrices in 3 hours with a sensitivity of 10 pg/mL/cytokine. Human and Mouse Cytokine Assay Kits are available in 1- and 10-plate sizes for fixed, off-the-shelf formats, or By Request, customer defined mix-to-order formats. By Request orders are processed and delivered in a pre-mixed ready to use format in ~1 week. Procarta Cytokine Assay Kits contain all the components to process cell culture supernatant samples and include: pre-mixed, ready to use, antibody-conjugated beads; 96-well filter plate and holder; assay and wash buffers; sample buffer; pre-mixed, ready to use detection antibodies; Streptavidin-PE (SAPE) fluorescence detection reagent; and premixed antigen standards

### How It Works


Procarta protein assays use the xMAP® technology (multi-analyte profiling beads) to enable the detection and quantitation of multiple protein targets simultaneously. The xMAP system combines a flow cytometer, fluorescent-dyed microspheres (beads), lasers and digital signal processing to effectively allow multiplexing of up to 100 unique assays within a signal sample. The Procarta assay kits are compatible with all Luminex and Luminex-based instruments currently available.




**1** Incubate your sample with the antibody-conjugated beads for 30 minutes




**2** Add detection antibody and incubate for 30 minutes



**3** Add SAPE and incubate for 30 minutes



**4** Detect interactions on a Luminex instrument



- Profile Cytokines in Diverse Matrices**

Quantitate cytokines from cell culture supernatants, serum or plasma samples

- Minimal Sample Requirements**

Results from up to 28 cytokines with 25  $\mu$ L of sample

- Complete Assay in 3 Hours**

Sample to quantitative results with 10 pg/mL/ cytokine sensitivity

- Flexible Kit Formats**

Pre-formatted or By Request Cytokine Assays delivered within ~1 week

### Specifications

Sensitivity	10 pg/mL/cytokine
Precision	Average Inter-assay CV <10% Average Intra-assay CV <10%
Spike Recovery (Accuracy)	80-120%
Cross Reactivity	Negligible
Matrices	Cell culture supernatants, serum, plasma



## Measure Protein and Gene Expression from the Same Sample Well


Panomics' Procarta Cytokine Assay Kits enables the profiling of up to 26 different cytokines per reaction. When coupled with Panomics' QuantiGene Plex assays, both gene and protein expression can be quantitated from the same sample well (supernatant for protein quantitation and cell lysate for mRNA quantitation) enabling the parallel study of gene expression at the RNA and protein level.

## Procarta Standard Diluent Kits

Procarta Cytokine Standard Diluent Kits for plasma and serum, sold separately, contain a single component, species/matrix-specific Standard Buffer, for preparation of antigen standards and dilution of experimental samples

(if required). Panomics' Procarta Cytokine Standard Diluent Kits are designed for use with Procarta Cytokine Assay Kits. Using the appropriate Standard Buffer ensures optimal recovery and sensitivity of the cytokines being analyzed in a serum or plasma matrix.

Human 28	Mouse 22	Rat 8
IL-1 $\beta$	IL-1 $\alpha$	IL-1 $\alpha$
IL-2	IL-1 $\beta$	IL-1 $\beta$
IL-4	IL-2	IL-6
IL-5	IL-3	TNF- $\alpha$
IL-6	IL-4	MCP-1
IL-7	IL-5	MIP-1 $\alpha$
IL-8	IL-6	KC
IL-10	IL-10	RANTES
IL-12(p70)	IL-12(p40)	
IL-12(p40)	IL-12(p70)	
IL-13	IL-17	
IL-17	IL-13	
IFN $\gamma$	KC	
GM-CSF	RANTES	
TNF- $\alpha$	IFN $\gamma$	
TNF- $\beta$	GM-CSF	
G-CSF	TNF- $\alpha$	
MIP-1 $\beta$	MIP-1 $\alpha$	
MIP-1 $\alpha$	EOTAXIN	
MCP-1	MCP-1	
EOTAXIN	MCP-3	
FGF basic	VEGF	
VEGF		
RANTES		
Leptin		
ENA-78		
NGF		
IP10		



Standard Curves were generated in cell culture media using Panomics' 19-plex Procarta Mouse Cytokine Assay Kit. Each analyte has a sensitivity of 10 pg/mL or less and an assay range over 3 logs.

## Procarta Cytokine Assay Kits, Preformatted

Product	Species	Plex	Size	Cat. #
Procarta Cytokine Assay Kit	Human	28	96 assays (1-plate)	PC0012
Procarta Cytokine Assay Kit	Human	20	96 assays (1-plate)	PC0001
Procarta Cytokine Assay Kit	Human	10	96 assays (1-plate)	PC0003
Procarta Cytokine Assay Kit	Mouse	19	96 assays (1-plate)	PC0101
Procarta Cytokine Assay Kit	Mouse	10	96 assays (1-plate)	PC0103
Procarta Cytokine Assay Kit	Mouse	10	96 assays (1-plate)	PC0103
Custom Blend Assay Kit		X-Plex	96 assays (1-plate)	Inquire

Check our website for additional sizes and custom blend options

## Procarta Cytokine Standard Diluent Kits

Product	Species	Matrix	Size	Cat. #
Procarta Cytokine Standard Diluent	Human	Serum	2.5 mL	PC0050
Procarta Cytokine Standard Diluent	Human	Plasma	2.5 mL	PC0052
Procarta Cytokine Standard Diluent	Mouse	Serum	2.5 mL	PC0150
Procarta Cytokine Standard Diluent	Mouse	Plasma	2.5 mL	PC0152

Check our website for additional sizes



### U.S. Corporate Headquarters

Panomics, Inc.  
6519 Dumbarton Circle  
Fremont, CA 94555  
Toll Free: 877 PANOMICS (1.877.726.6642)  
Direct: 1.510.818.2600  
Fax: 1.510.818.2610

Email: [info@panomics.com](mailto:info@panomics.com)

Email: [orders@panomics.com](mailto:orders@panomics.com)

Email: [techsupport@panomics.com](mailto:techsupport@panomics.com)

### European Headquarters

Panomics Srl  
Via Sardegna 1  
20060 Vignate-Milano (Italy)